

# DISPOSABLE ELECTRODES FROM A NAME YOU TRUST!



All Aaron® blade, ball, needle, loop, and square disposable electrodes use the standard 3/32" diameter stainless steel shafts. These electrodes, like all Aaron® electrodes, feature safety grip insulators, combining patient and user safety with easy insertion into and removal from the surgical pencil.

The loop and square electrodes feature tungsten wire for superior shape and wire integrity throughout the excision procedure, while minimizing thermal damage and yielding excellent histopathology specimens. Each molded loop and square electrode includes our patented tray to protect the loops from bending during shipping and handling.

Ball electrodes in 3mm and 5mm sizes are available for fulguration and desiccation, while needle electrodes are available for precise, pinpoint cutting and coagulation.

SuperCut™ Needles (ES 60-63) are made of Tungsten and keep their shape throughout the procedure.

The Aaron Arthroscopic electrodes are available in four tip styles and designed to fit all standard electro-surgical pencils.



7100 30th Avenue North  
 St. Petersburg, FL 33710-2902  
 U.S. Ph 1-800-537-2790 • Fx 1-800-323-1640  
 Int'l. Ph +1-727-384-2323 • Fx +1-727-344-3876  
 www.aaronmed.com • sales@aaronmed.com

# LIST OF STYLES AND CONVERSION CHART

		STRYKER/COLORADO	PREMIER	WALLACH	CONMED	VALLEY	ACMI/CHICON	UTAH
SuperCut™ Tungsten Needles	ES60	3CM Super Fine Needle, 5/box	N101A					UTAH
	ES61	3CM Modified Super Fine Needle, 5/box	N103A			E1651		
	ES62	4.5CM Super Fine Needle, 5/box	N102A					
	ES63	4.5CM Modified Super Fine Needle, 5/box	N104A					
Needle Electrodes	ES02	Standard Needle, 25/box	COOPER	PREMIER	WALLACH	CONMED	VALLEY	UTAH
	ES03	Extended Needle, 25/box			909135	275-1502/9113	E1552	
	ES38	Modified Needle, 25/box				275-1504/9116	E1552-6	
	ES40	Extended Modified Needle, 25/box				138025		6162901
Dermal Loops	ES22	1/4" Short Shaft Loop, 5/box						
	ES23	5/16" Short Shaft Loop, 5/box						
	ES24	3/8" Short Shaft Loop, 5/box						
	ES25	1/2" Short Shaft Loop, 5/box						
	ES26	3/4" Short Shaft Loop, 5/box						
	ES26	3/4" Short Shaft Loop, 5/box						
Dermal Tips	A804	Sharp Dermal Tip, Non-Sterile, 100/box						
	A806	Blunt Dermal Tip, Non-Sterile, 100/box						
Blade Electrodes	A807	Blunt Dermal Tip, Sterile, 50/box						
	ES01	Standard Blade, 50/box				275-1500/9111	E1551G	
	ES04	Extended Blade, 25/box				275-1503/9114	E1551-6	
	ES37	Modified Blade, 25/box				275-1524	E1553	
	ES39	Extended Modified Blade 25/box				275-1525		
	ES18	Anoled Blade, 5/box						
Ball Electrodes	ES06	Extended 3mm Ball, 5/box	B0312					DBL-311
	ES07	Extended 5mm Ball, 5/box	B0512			60-5181-102	E1563	6156901
	ES20	3/16" Ball, 25/box				60-5181-103	E1564	6157901
	ES21	5/32" Ball, 25/box				275-1501/9112	E1550	DBL-511
	ES50	1/8" Ball, 25/box						
LLETZ Loop Electrodes	ES08	5mm x 5mm Loop, 5/box						DL-P-T11
	ES09	10mm x 10mm Loop, 5/box	R1010					DL-P-S11
	ES10	15mm x 8mm Loop, 5/box	R1507*					
	ES11	15mm x 10mm Loop, 5/box						
	ES12	20mm x 8mm Loop, 5/box	R2008			60-5180-102		DL-P-M11
	ES13	20mm x 15mm Loop, 5/box	R2015					DL-P-L11
	ES31	20mm x 20mm Loop, 5/box	R2020					
	ES42	20mm x 12mm Loop, 5/box						DL-P-W11
	ES43	15mm x 15mm Loop, 5/box						
	ES44	15mm x 5mm Loop, 5/box	R1515					
	ES44	15mm x 5mm Loop, 5/box	R1505					
	ES45	13mm x 13mm Loop, 5/box						
	ES46	10mm x 5mm Loop, 5/box						
	ES47	25mm x 10mm Loop, 5/box						
	ES48	10mm x 3mm Loop, 5/box	R2510					
	ES49	20mm x 10mm Loop, 5/box						
ES51	13mm x 8mm Loop, 5/box	R2010			60-5180-104			
ES52	10mm x 8mm Loop, 5/box				60-5180-101*			
ES53	15mm x 12mm Loop, 5/box							
ES57	25mm x 8mm Loop, 5/box							
ES58	25mm x 22mm Loop, 5/box							
LLETZ Square Electrodes	ES14	5mm x 5 mm Square, 5/box						E1565
	ES15	10mm x 4mm Square, 5/box	S1004					E1566
	ES16	10mm x 8mm Square, 5/box						E1567
	ES41	5mm x 10mm Square, 5/box						6155901
	ES56	10mm x 10mm Square, 5/box						
Arthroscopic Electrodes	AR00	Electrode Hook 90, 5/box						E1510
	AR01	Electrode Hook 45, 5/box						E1511
	AR03	Electrode Meniscotomy, 5/box						E1512
								E1513

\* shows 1mm difference in size

## Loop Size Locator

